









# Peel Away Sheath Introducer



CODE	SIZE	COLOUR	DILATOR LENGTH (MM)	DILATOR EXTERNAL DIAMETER (MM)	SHEATH LENGTH (MM)	SHEATH INTERNAL DIAMETER (MM)
432005	5 F		180	1.67	142	1.86
432006	6 F		180	2.05	142	2.17
432007	7 F		180	2.37	142	2.51
432008	8 F		180	2.67	142	2.84
432009	9 F		180	3.05	142	3.16
432010	10 F		180	3.35	142	3.48
432011	11 F		180	3.66	142	3.80
432012	12 F		180	4.00	142	4.14
432013	13 F		180	4.33	142	4.50
432014	14 F		180	4.66	142	4.80
432016	16 F		180	5.40	142	5.60
432018	18 F		180	6.00	142	6.20
432020	20 F	180	6.70	142	6.90	



ISOMED Peel Away Introducers allow for smooth catheter insertion and manipulation. The X.R.O. peel away sheath allows the device to be completely removed without shifting the catheter out of position.

The dilator has a long rounded tip for atraumatic insertion. The taper from the dilator to the sheath is precise and smooth, allowing easy sheath introduction.

Supplied sterile (ethylene oxide), with a J guidewire of 0.035" x 45 cm and a puncture needle of 18G x 75 mm.

Individual puncture needles, dilators and sheaths in different diameters and lengths are available upon request.