










## Biliary Catheters - Single & Double Lumen



CODE	COLOUR	EXTERNAL DIAMETER	LENGTH	MAX. INFLATE	BALLOON DIAMETER (MAX.)	
<b>BILIARY CATHETERS SINGLE LUMEN</b>						
340523		White	5F (1,70 mm)	23 cm	1,50 ml	10 mm
340540		White	5F (1,70 mm)	40 cm	1,50 ml	10 mm
340623		Blue	6F (2,00 mm)	23 cm	2,00 ml	13 mm
340640		Blue	6F (2,00 mm)	40 cm	2,00 ml	13 mm
<b>BILIARY CATHETERS DOUBLE LUMEN</b>						
540540		White	5F (1,70 mm)	40 cm	1,50 ml	10 mm
540740		White	7F (2,35 mm)	40 cm	2,50 ml	14 mm

Other lengths available upon request



ISOMED Biliary Catheters are made with X.R.O. polyurethane. They are colour coded for easy identification. The specially designed latex balloon, fixed with two sutures, is perfectly symmetrical. Biliary catheters are graduated every centimeter

from the distal end, and are re-inforced with a removable steel wire. The double lumen catheters allow the injection and opacification during the procedure.

Double packaged - Sterile (ethylene oxide)