




## Huber Needles



CODE	COLOUR	DIAMETER		LENGTH
		MM	GAUGE	MM
HD 2220	Black	0,7	22	20
HD 2225				25
HD 2230				30
HD 2238				38
HD 2260				60
HD 2020	Yellow	0,9	20	20
HD 2025				25
HD 1920	Off white	1,1	19	20
HD 1925				25



Packaging : 50 units box - Sterile Ethylene Oxide



**FB Medical Huber Needles** are available with a luer lock hub. Their unique point design promotes easy penetration and maximum life for the port septum. They are non-siliconized to minimise "back out".

CODE	COLOUR	DIAMETER		LENGTH
		MM	GAUGE	MM
HC 2220	Black	0,7	22	20
HC 2225				25
HC 2230				30
HC 2020	Yellow	0,9	20	20
HC 2025				25
HC 2030				30
HC 1920	Off white	1,1	19	20
HC 1925				25
HC 1932				32



Packaging : 50 units box - Sterile Ethylene Oxide



**FB Medical** offers a complete range of **Huber Needles, straight or 90° curved**. They are designed for single use and easy insertion into vascular access device, and may be used for continuous or bolus injections of pain control drugs, chemotherapeutics agents, parenteral nutrition and blood products.