









Huber Needles



CODE	COLOUR	Ø		LENGTH		
		GAUGE	MM	MM	INCH	
HD2225		Black	22	0.7	25	1"
HD2238		Black	22	0.7	38	1 1/2"
HD2260		Black	22	0.7	60	2 3/8"
HD2025		Yellow	20	0.9	25	1"
HD1925		Off White	19	1.1	25	1"



FB medical Huber Needles are available with a luer lock hub. Their unique point design promotes easy penetration and maximum life for the port septum. They are non-siliconized to minimise "back out".

CODE	COLOUR	Ø		LENGTH		
		GAUGE	MM	MM	INCH	
HC2220		Black	22	0.7	20	3/4"
HC2225		Black	22	0.7	25	1"
HC2230		Black	22	0.7	30	1 1/4"
HC2025		Yellow	20	0.9	25	1"
HC2030		Yellow	20	0.9	30	1 1/4"
HC1925		Off White	19	1.1	25	1"

FB Medical offers a complete range of Huber needles, straight or 90° curved. They are designed for single use and easy insertion into vascular access device, and may be used for continuous or bolus injections of pain control drugs,

chemotherapeutics agents, parenteral nutrition and blood products.

Packaging : 50 units box - Sterile Ethylene Oxide